

# SUPERNOVA Stent Retriever



### Advanced Clot Capture Technology



### Spiral Fusiform Technology

Ensures superior engagement and retrieval of both red and hard fibrin-rich (white) clots.



### Optimized Radial Force

Provides effective deployment and retrieval, maintaining consistent performance.



### Radiopaque Markers

Ensures visual confirmation under fluoroscopy for accurate positioning.



### Dynamic Strut Activity

Engineered for optimal flexibility and control during deployment and retrieval.

The Gravity Supernova Revascularization Device is a stent-based thrombus retrieval device with a collapsible, fully retrievable, self-expanding laser cut stent which is mounted on a wire shaft for delivery.

## Designed for Precision and Performance



Achieved **95% mTICI  $\geq$  2B recanalization success rate**, with 64% first pass success in first-in-human clinical studies



Innovative design enhances the capture and removal of **hard fibrin-rich (white) clots**

### Supernova Revascularization Device Information

Model	Working Length	Stent Diameter	Vessel Diameter Range	Radiopaque Markers	Delivery Catheter ID
GRVY-004-043	39mm	4mm	1.5mm - 4.0mm	12	0.021"
GRVY-006-047	43mm	6mm	2.0mm - 5.5mm	15	0.021"



### Optimize Your Stroke Intervention

Join the leading experts in stroke care and enhance your clinical practice with the Supernova Stent Retriever. Experience the unmatched performance and reliability that Supernova offers in the most challenging clinical scenarios.

**CAUTION:** These devices are restricted to sale, distribution, and use by or on the order of a physician. Detailed indications, contraindications, warnings, and instructions for use are provided in the accompanying Instructions for Use (IFU).

**INDICATIONS FOR USE:** The Supernova Revascularization Device is indicated for temporary use to restore blood flow in the cerebral vasculature of patients suffering from an acute ischemic stroke.

**HOW SUPPLIED:** The Supernova Revascularization Device is provided sterile (Ethylene Oxide, non-pyrogenic), as a single unit, and should be stored in a dry, cool place.

**DISCLAIMER:** The content provided is for informational purposes only and should not be construed as medical advice. Users should consult with a qualified healthcare professional before making any medical decisions based on this information. Gravity Medical Technology, Inc. assumes no responsibility for any errors or omissions in the content.

© 2024 Gravity Medical Technology, Inc. All Rights Reserved.



105 N. Pointe Drive, Suite E,  
Lake Forest, CA 92630  
United States of America

+1 (949) 328-6609

admin@gravitymedtech.com

www.gravitymedtech.com